RANGER[®] 7000

RS232 Interface
Instruction Manual

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
INTRODUCTION

This RS232 interface Kit is for use with the Ohaus Ranger[®] 7000.

KIT CONTENTS

- Pressure Button
- Screws (4)
- RS232 PC Board


Note: RS232 cable is not included. If RS232 cable is needed from OHAUS, please order 80500552.

INTERFACE INSTALLATION

Only one interface module can be installed in the terminal option slot at a given time. Install the Interface module on the underside of the terminal as shown below.

Note: The appearance of your model may be different from the one shown.

1. Release the two screws, remove the red interface options cover, and unplug the black rubber insert.

-
- Diagram of the back of the device showing the internal compartment. Labels include: Cord groove, Screws, Pins.

- ## SETUP

- A) The OHAUS Interface is preset to communicate using the following settings: 9600 baud, 8 bit, no parity, no handshake. Determine the RS232 parameters required for the printer or computer that is to be connected. See computer or printer documentation if assistance is required. If the parameters do not match, it will be necessary to change either the scale's settings, or the computer / printer settings.
- B) Upon installation of the interface, the scale will recognize the RS232 Interface and add relevant items to the menu. Configure the scale to the desired RS232 and printing parameters; refer to the scale's Instruction Manual for assistance in using the menus.

Communication

RS232

USB

2nd RS232

▼

▲

↻

✓

Configuration

Reset

Baud Rate

9600

Parity

8 None

Stop Bits

1 bit

Handshake

None

ASCII Alternate Command

▼

▲

↻

✓

Print Setup

Reset

Stable Weight Only

Off

Auto Print

Print Template

Simple

Edit Template

Simple

Line Feed

1 Line

▼

▲


↻

✓

CONNECTION

The RS232 Interface Kit cable terminates with a 3 pin Sub-D female connector.

Active pins: Pin 1 = Ground, Pin 2 = RXD, Pin 3 = TXD

The Interface connector will connect directly to a computer, an Ohaus printer and many other printers

Note: When the second RS232 interface has been installed, an icon will appear on the display:

OPERATION

RS232 OUTPUT

With the interface installed the terminal will operate in one of the following two ways:

- Auto Print = off, Stable = on or off
Pressing PRINT will send the display data.
If Stable is set to On, the balance will wait for a stable reading before sending the stable data.
- Auto Print = not off, Stable = on or off
The balance will automatically send data based on the settings in the Menu.
If Stable is set to On, only stable values will be sent.

Data sent from the Interface is in standard ASCII format terminated with a [CRLF] (carriage return-line feed). The output format is as follows:

✓ Printout string for g, kg, lb, oz units:

[weight]	9 characters (right justified)
[space]	1 character
[unit]	3 characters max (left justified)
[space]	1 character
[stability indicator]	1 character "?" when unstable, blank when stable
[space]	1 character
[G/N]	1 character
[space]	1 character
[Miscellaneous]	5 character
[Term. Char(s)]	2 character

✓ Printout string print format for lb:oz unit

[weight]	4 characters (right justified)
[space]	1 character
[unit]	2 characters max (left justified)
[space]	1 character
[weight]	7 characters (right justified)
[space]	1 character
[unit]	2 characters max (left justified)
[space]	1 character
[stability indicator]	1 character "?" when unstable, blank when stable
[space]	1 character
[G/N]	1 character
[space]	1 character
[Miscellaneous]	5 character
[Term. Char(s)]	2 character


Output example:

**Weigh Mode
Printout**

	Description	Comment
11.11 kg NET	Result line	If Printx → Content → Result is ON
12.34 kg G	Gross value line	If Printx → Content → Gross is ON and a tare value is entered
11.11 kg NET	Net value line	If Printx → Content → Net is ON and a tare value is entered
1.23 kg T	Tare value line	If Printx → Content → Tare is ON and a tare value is entered
MODE: WEIGH	Mode line	If Printx → Content → Application Mode is ON
<no line printed>	Information line	If Printx → Content → Into is ON

COMPLIANCE

Disposal

In conformance with the European Directive 2002/96 EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements. Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.

Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

For disposal instructions in Europe, refer to www.ohaus.com/weee.

Thank you for your contribution to environmental protection.

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The image shows two side-by-side screenshots of a scale's configuration menu. The left screen is titled 'Kommunikation' (Communication) and lists 'RS232', 'USB', and 'RS-232 optional' with right-pointing arrows. 'RS-232 optional' is highlighted. The right screen is titled 'Konfiguration' (Configuration) and lists 'Zurücksetzen' (Reset), 'Baudrate' (9600), 'Parität' (None), 'Stoppbits' (1 bit), 'Handshake' (None), and 'ASCII Alternativer Befehl' (ASCII Alternative Command), all with right-pointing arrows. 'Zurücksetzen' is highlighted. At the bottom of both screens are four navigation buttons: down arrow, up arrow, left arrow, and a checkmark.

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