

INSTRUCTIONS

Eye Wash Bottle Wall Rack

Catalog No. 24859-0000



The epoxy coated steel rack holds a single 32 ounce (1000cc) Eye Wash Bottle, F24851-0000. It is designed to be hung on a wall by means of a single screw or nail. It is important that the screw be anchored properly, into a wooden stud or into a plastic or lead anchor if no stud is available. If a nail is used it must go into a solid wooden member. Other anchoring systems may be used for wallboard or plaster walls.

OTHER FINE SCIENCEWARE® PRODUCTS FOR YOUR LABORATORY

Eye Wash Station

Catalog No. 24865-0000, 24866-0000



Bright, Easy to Locate

This easily recognized, ready-to-use Eye Wash Station comes complete with eye wash bottle, dust cover, mirror, and instructions for use and service. A bracket holds the bottle securely in place, yet the bottle is easily released when pulled out. Backboard is 44 x 38cm (17½" x 15") with 4 screw holes for wall mounting. Includes station and bottle only, eye wash not included. Bottles supplied empty.

Catalog No.
F24865-0000
F24866-0000

Bottle Size
Station with 500ml (16oz) empty eye wash bottle
Station with 1000ml (32oz) empty eye wash bottle